

# Square Storage Bottles & Carboys - PETG

## PETG Storage Bottles



**PETG Square Storage Bottles – are Durable, Ideal for Storage and Sampling of Culture Media, General Laboratory Reagents and Biological Buffers.**

PETG's enhanced material properties features low permeability to  $CO_2/O_2$  & provide longer term storage of sensitive contents



### Details:

- Made from durable, break resistant PETG suitable for frozen storage down to -40C
- Leakproof bottles and closures are provided sterile and non pyrogenic, eliminating washing, depyrogenation steps.
- Comprised of resins that are ADCF free
- Supported by Validation documents
- Packed in shrink wrapped trays

Cat No.	Product Description	Volume	Pack Size
GX-590005	PETG Bottle Steile, Clear bottle, white cap,	5ml	500/CS
GX-590010	PETG Bottle Steile, Clear bottle, white cap,	10ml	500/CS
GX-590030	PETG Bottle Steile, Clear bottle, white cap,	30ml	280/CS
GX-590060	PETG Bottle Steile, Clear bottle, white cap,	60ml	200/CS
GX-590125	PETG Bottle Steile, Clear bottle, white cap,	125ml	96/CS
GX-590250	PETG Bottle Steile, Clear bottle, white cap,	250ml	60/CS
GX-590500	PETG Bottle Steile, Clear bottle, white cap,	500ml	40/CS
GX-591000	PETG Bottle Steile, Clear bottle, white cap,	1000ml	24/CS

# Square Storage Bottles & Carboys - PETG

## PETG Carboys



Carboy is an advanced fluid handling system developed with the end-user in mind. The ergonomically designed body is coupled with FlexiCap technology to provide the user an easy-to-use and easy-to-clean carboy. Carboy is made from a clear, BPA-free, USP Class VI, FDA Grade resin. This Carboy is suitable for aqueous solutions, as well as most life sciences buffer solutions.

Also Available in sterile format



2.5L (0.5 Gal) PETG Carboy with 83B Cap



5L (1 Gal) PETG Carboy with 83B Cap



10L (2.5 Gal) PETG Carboy with 83B Cap



20L (5 Gal) PETG Carboy with 120mm Cap

### Details

- Novel shape maximizes bench and shelf space
- Bottom handles provide a secure grip while lifting and transporting
- Large top handles for a full hand grip
- Molded carboy identification marks
- Large, clear graduation marks with  $\pm 5\%$  accuracy
- FlexiCap™ technology with versatile adapters and replaceable gasket
- EZBoy™ body design allows for maximum solution recovery
- Available in 2.5L, 5L, 10L, 20L capacities

### Caution

- Do not use carboys under pressure or full vacuum
- Do not use carboys at temperatures greater than the rated maximum
- Do not store oxidizing agents in carboys
- Do not flame carboys
- Do not use carboy spigot as a handle
- Do not expose to organic compounds.
- Carboy will become prone to stress fracturing

Cat No.	Product Description	Volume	Pack Size
GX-851-0111	Carboy, PETG FlexiCap™ 80	2.5L	1/Pk
GX-851-1111	Carboy, PETG FlexiCap™ 80	5L	1/Pk
GX-851-2111	Carboy, PETG FlexiCap™ 80	10L	1/Pk
GX-851-3121	Carboy, PETG FlexiCap™ EZFill™ 120	20L	1/Pk



Genetix Biotech Asia Pvt. Ltd.

71/1, First Floor, Shivaji Marg, Najafgarh Road, New Delhi - 110015 - INDIA  
 phone : +91-11-45027000 ~ email : info@genetixbiotech.com ~ www.genetixbiotech.com