



- Round Carboys are manufactured using high purity virgin USP Class VI Medical grade polypropylene
- These carboys are used for storing and shipping reagents, chemicals, culture media, and standards
- PP Carboys are autoclavable at 121°C for twenty minutes
- Can be used at a temperature below 0° C
- Superior chemical resistance

Caution

- Autoclave carboy, cap, and spigot separately

A premium line of high quality, round polypropylene (PP) laboratory carboys with the FlexiCap® technology. These durable, rugged and shatterproof round carboys are a great alternative to glass carboys. They can be used for mixing, storing, and dispensing liquids in volumes from ten to twenty liters. These carboys are available with and without a spigot.

Round Carboys are compatible with FlexiCap® that is leak-proof, has replaceable platinum-cured silicone gaskets, and provides the ultimate flexibility. The open FlexiCap® allows a variety of configurations by changing the adapter insert, enabling users to quickly change the top connections of each carboy or bottle to match their application needs.



4

5



1 83mm FlexiCap® Technology

2 Handle design to make handling carboy easier

3 Round carboy design maximizes mixing efficiency

4 Easy to read, graduation marks on the twenty liter

5 Spigot allows simple and accurate dispensing of media and solutions



Product Number	Description	Capacity (L)	Cap Size (mm)	Color	Autoclavable	Qty (EA)
GX-551-711B	10L Round PP Carboy with 83B Cap	10	83	Opaque	Yes	1
GX-551-711C	20L Round PP Carboy with 83B Cap	20	83	Opaque	Yes	1
GX-551-718B	10L Round PP Carboy with 83B Cap and Spigot	10	83	Opaque	Yes	1
GX-551-718C	20L Round PP Carboy with 83B Cap and Spigot	20	83	Opaque	Yes	1



Laboratory Carboys

We offer a premium line of high quality, innovative laboratory carboys. Our versatile carboys aim to take full advantage of storage and/or transportation efficiency at minimal effort. The shape, size of the handle, grip, and neck, and imprints were designed and developed to provide the ultimate end-user experience. Furthermore, all laboratory carboys and bioprocess carboys are BPA-free, made from a USP Class VI material, and are FDA Grade for optimum protection of your assets. Materials available include; HDPE, PP, PETG, PC, Amber HDPE, Dark Amber PP, and options with or without spigot and sanitary top!