

DATA SHEET

PHOSPHATE BUFFER SALINE (PBS) TABLETS

Cat# PG-404-100TABS

Phosphate Buffered Saline Tablets (PBS) are a balanced salt mixture for easy PBS buffer preparation for a variety of applications including reagent preparation, diluting cells for flow cytometry and cell culture reagent.

No need for weighing and mixing chemicals

Simply dissolve PBS tablets into water

INSTRUCTION OF USE

Dissolve each tablet of PBS with 100 ml of distilled water. This buffer contains 10 mM phosphate, 150 mM sodium chloride, pH 7.3 to 7.5.

STORAGE: Store at room temperature. Keep the tablets in a dry and cool place.

REMARKS: No preservative is added.