

## Agar Bacteriological, 500g

Cat No. PG-736-500G

CAS : 9002-18-0

Test	Specifications
Appearance (Colour)	Off - white to pale yellow
Appearance (Form)	Powder
Solubility (Turbidity) 1.5% aq. hot solution	Almost clear to slight hazy
Loss on drying	max. 15%
Sulphated Ash	max. 5%
Gel strength (1.5% gel)	min. 600 gm/cm <sup>2</sup>
Foreign organic matter	max. 1%
Starch	Passes test
Gelatin	Passes test
Water absorption capacity	Passes test