

## DATA SHEET

---

### **PHOSPHATE BUFFER SALINE (PBS) TABLETS**

**Cat# PG-404-100TABS**

Phosphate Buffered Saline Tablets (PBS) are a balanced salt mixture for easy PBS buffer preparation for a variety of applications including reagent preparation, diluting cells for flow cytometry and cell culture reagent.

No need for weighing and mixing chemicals

Simply dissolve PBS tablets into water

#### **INSTRUCTION OF USE**

Dissolve each tablet of PBS with 100 ml of distilled water. This buffer contains 10 mM phosphate, 150 mM sodium chloride, pH 7.3 to 7.5.

**STORAGE:** Store at room temperature. Keep the tablets in a dry and cool place.

**REMARKS:** No preservative is added.