

10X T4 DNA Ligase Buffer

Cat. No.

PGB69

Pack Size

1.5ml

Description

10X T4 DNA Ligase buffer is used with T4 DNA Ligase.

10X Buffer Composition

400 mM Tris-HCl, 100 mM MgCl₂, 100 mM DTT, 5 mM ATP
(pH 7.8 at 25°C)

CERTIFICATE OF ANALYSIS

Tested for the absence of endo-, exodeoxyribonucleases,
ribonucleases and functionally tested in DNA ligation.

