



PUREGENE

Water, Sterile, Nuclease-Free

PRODUCT CODE: PG-674-500ML

Puregene's sterile, Nuclease-free water is prepared using diethylpyrocarbonate (DEPC) to ensure the removal of RNase contamination. Ideal for molecular biology applications.

- RNase-free
- Sterile
- Optimized for molecular biology applications

Grade: BIOTECHNOLOGY GRADE

CAS Number: 7732-18-5

Synonym(s): None

Molecular Weight: 18.02

Molecular Formula: H₂O

Storage Condition: ROOM TEMPERATURE

Hazard Statements: NA

Precautionary Statements: NA

Specification:

DNase (P/F)NONE

Protease (P/F).....NONE

RNase (P/F)NONE

Sterility (P/F)..... PASS