

Oligo(dT)₁₈

Store at -20°C

Cat. No.

PGPO13

Pack Size

60ul (30ug)

5'-d (TTTTTTTTTTTTTTTTTTTT)-3'

Concentration: 100uM (0.5ug/ul)

Description

The oligo (dT)₁₈ primer is a synthetic single-stranded 18-mer oligonucleotide with 5'- and 3'-hydroxyl ends. The primer is supplied as a ready-to-use, 20X concentrated aqueous solution.

Application

First strand cDNA synthesis.

CERTIFICATE OF ANALYSIS

Functionally tested in RT-PCR using 100 fg of control GAPDH RNA.

