

250cards on 300gsm paper



FEATURES:


- **Broad range** – 10-200 kDa.
- **Ready-to-use** – supplied in a loading buffer for direct loading on gels.
- **Sharp bands**, reference band included.
- **Includes a 50 kDa reference band** of a greater intensity.

APPLICATIONS:

- Accurate protein sizing on SDS-polyacrylamide gels and Western blots.

Storage
Stable at +4°C for 3 months
Store at -20°C for 12 months

QUALITY CONTROL:
Tested in SDS-polyacrylamide gel electrophoresis and Western blotting.



8-16 % Tris-glycine SDS-PAGE stained with PageBlue Protein Staining Solution.
Recommended product amount per well in minigel applications is 7µl.



MFG
71/1, First Floor, Shivaji Marg, New Delhi - 110015 INDIA
Phones : +91-11-45027000, Fax : +91-11-2541-9631
E-mail : info@genetixbiotech.com Website: www.genetixbiotech.com



UNSTAINED PROTEIN LADDER

Cat # PG-PMT0772 **Size: 2X250µl**

Store at -20°C