

ERLENMEYER FLASKS - PC & PETG

Erlenmeyer Flasks are suitable for the culture of cell lines with high requirements on oxygen. It can also be used for the culture of bacteria, fungi and animal and plant cells in suspension.



Baffled Base



Flat Base

Product Description

- Two different cap types : Cover and breathable
- Made from optically clear polycarbonate or PETG
- Plain bottom and baffled bottom, baffled bottom improve mixing and increase gas exchange.
- Available with 4 volumes of 125ml, 250ml, 500ml and 1000ml
- Ideal for shaker culture applications and storage
- Molded-in graduations for accuracy
- Vent cap option for continuous gas exchange while ensuring sterility and preventing leakage
- Sterilized by gamma irradiation
- RNase-/DNase-free , Non-pyrogenic
- Temperature range for PC Flasks: -80°C to +121°C
- Temperature range for PETG Flasks: -60°C to +60°C



Ordering Information

Cat No	Items	Material	Cap Type	Application	Volume ml	Pack/Pc	Case/Pc
GX-PG125SS	Erlenmeyer Flasks / Shake Flasks plain bottom	PETG	seal caps	cell culturing/ bacteria culturing	125	1	50
GX-PG125FS	Erlenmeyer Flasks / Shake Flasks plain bottom	PETG	filter caps	cell culturing/ bacteria culturing	125	1	50
GX-PG250SS	Erlenmeyer Flasks / Shake Flasks plain bottom	PETG	seal caps	cell culturing/ bacteria culturing	250	1	50
GX-PG250FS	Erlenmeyer Flasks / Shake Flasks plain bottom	PETG	filter caps	cell culturing/ bacteria culturing	250	1	50
GX-PG500SS	Erlenmeyer Flasks / Shake Flasks plain bottom	PETG	seal caps	cell culturing/ bacteria culturing	500	1	25
GX-PG500FS	Erlenmeyer Flasks / Shake Flasks plain bottom	PETG	filter caps	cell culturing/ bacteria culturing	500	1	25
GX-PG1000SS	Erlenmeyer Flasks / Shake Flasks plain bottom	PETG	seal caps	cell culturing/ bacteria culturing	1000	1	25
GX-PG1000FS	Erlenmeyer Flasks / Shake Flasks plain bottom	PETG	filter caps	cell culturing/ bacteria culturing	1000	1	25
GX-PC125SS	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	seal caps	cell culturing/ bacteria culturing	125	1	50
GX-PC125FS	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	filter caps	cell culturing/ bacteria culturing	125	1	50
GX-PC250SS	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	seal caps	cell culturing/ bacteria culturing	250	1	50
GX-PC250FS	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	filter caps	cell culturing/ bacteria culturing	250	1	50
GX-PC500SS	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	seal caps	cell culturing/ bacteria culturing	500	1	25
GX-PC500FS	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	filter caps	cell culturing/ bacteria culturing	500	1	25
GX-PC1000SS	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	seal caps	cell culturing/ bacteria culturing	1000	1	25
GX-PC1000FS	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	filter caps	cell culturing/ bacteria culturing	1000	1	25
GX-PC3000FS-AM	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	Vent caps	cell culturing/ bacteria culturing	3000	1	4
GX-PC5000FS-AM	Erlenmeyer Flasks / Shake Flasks plain bottom	PC	filter caps	cell culturing/ bacteria culturing	5000	1	4

Cat No	Item	Material	Cap Type	Volume ml	Pack Size
GX-011125B	Erlenmeyer Flasks / Shake Flasks Baffled Bottom	PC	seal caps	125	24
GX-012125B	Erlenmeyer Flasks / Shake Flasks Baffled Bottom	PC	filter caps	125	24
GX-011250B	Erlenmeyer Flasks / Shake Flasks Baffled Bottom	PC	seal caps	250	12
GX-012250B	Erlenmeyer Flasks / Shake Flasks Baffled Bottom	PC	filter caps	250	12
GX-011500B	Erlenmeyer Flasks / Shake Flasks Baffled Bottom	PC	seal caps	500	12
GX-012500B	Erlenmeyer Flasks / Shake Flasks Baffled Bottom	PC	filter caps	500	12
GX-011000B	Erlenmeyer Flasks / Shake Flasks Baffled Bottom	PC	seal caps	1000	24
GX-012000B	Erlenmeyer Flasks / Shake Flasks Baffled Bottom	PC	filter caps	1000	24



Genetix Biotech Asia Pvt. Ltd.

71/1, First Floor, Shivaji Marg, Najafgarh Road, New Delhi - 110015 INDIA
 Phone : +91-11-45027000 ~ Fax : +91-11-25419631 / E-mail : info@genetixbiotech.com
 For further detail please visit : www.genetixbiotech.com